



**DSS Studio**

---

## A Novel Mixed Reality Navigation System for Precise Tunnel Localization in ACL Joint Bone Puncture Surgery: A Navigation Plugin

Diyu Zhou ([dz372@cornell.edu](mailto:dz372@cornell.edu))

Shuhan Liang ([sl3537@cornell.edu](mailto:sl3537@cornell.edu))

Sissi Zheng ([wz472@cornell.edu](mailto:wz472@cornell.edu))

---

### Plugin Overview

This plugin, based on the Fologram framework, is designed as an AR visualization tool for bone marrow puncture guidance using inputs from different skeletal models. By optimizing the model and converting point clouds, it generates features such as bone marrow pathways, path guidance, distance and angle indicators, and dangerous area warnings. The goal is to enhance surgical precision and streamline the surgical process for orthopedic surgeons.

#### Supported Features (Version 1.0)

- Precise multi-target bone tunnel puncture point surgical navigation
- Real-time multi-target precise distance measurement
- Optimal path navigation
- Puncture angle calculation
- Implant distance and depth calculation
- Danger zone calculation

---

### Installation

bone\_tunnel

---

#### Software required

- Rhino 8.0
- Grasshopper
- Fologram

#### Hardware required

- PC
- HoloLens

#### Installation Step-by-step

- Download and open rhino8
- Download and open grasshopper
- Set up "bone\_tunnel" plug
- Set up "Fologram" plug


Guide of  
Functionality

Input




Spatial points  
3D skeletal model  
Bone spindle

Output





Surgical path  
Distance between bone tunnel points  
Bone tunnel depth  
Bone marrow canal opening angle

Feature Pseudo Code



- <sup>1</sup> Real-time multi-target optimal distance measurement & navigation
- <sup>2</sup> Puncture angle calculation
- <sup>3</sup> Implant depth calculation
- <sup>4</sup> Danger zone calculation

## Detail of Functionality



### Remesh Bone

This component is used to optimize the input skeletal mesh to facilitate subsequent calculations and optimizations.

#### **Input- edge\_length**

Type: float

##### Description

The `edge_length` parameter determines the length of the edges of the target mesh (measured in the current model units of Rhino).

This parameter controls the resolution of the mesh, affecting its density and level of detail.

##### Functionality:

A smaller `edge_length` generates a higher-resolution mesh (with more faces).

A larger `edge_length` generates a lower-resolution mesh (with fewer faces).

#### **Input-Preserve**

Type: bool (Boolean)

##### Description:

The `preserve` parameter determines whether to maintain the boundary shape of the input mesh (usually referring to sharp edges or boundary lines of the mesh).

##### Functionality:

When set to `True`, QuadRemesh will attempt to preserve the original mesh's boundary shape as much as possible.

When set to `False`, QuadRemesh will not enforce boundary preservation, allowing the mesh to be freely adjusted based on the target edge length.

#### **Input-mesh\_input**

Type: Rhino.Geometry.Mesh

##### Description:



The `mesh_input` parameter represents the original mesh to be remeshed, defining the geometry to be reconstructed.

### Output-mesh\_output

Type: Rhino.Geometry.Mesh

#### Description:

The output is an optimized mesh, making it convenient for subsequent point set conversion.



### Bone to Point

This component is used to convert the input skeletal mesh into a point cloud to facilitate subsequent calculations and optimizations.

### input-Bone

Type: Rhino.Geometry.Mesh

#### Description:



The input value can be context-dependent, either the original mesh bone or the optimized skeleton mesh.

### Output-Vertices

Type: Point

#### Description:

The output consists of vertices converted from the mesh skeleton.



### Tunnel Point Select

The selection of specific bone tunnel points on the skeletal mesh is done through the

point set on the skeletal mesh.

**input-vertices**

Type: Point 3D

Description:

This input value is the point set generated from the skeletal mesh.

**input-point\_a\_index**

Type: Int

Description:

Points are selected using the index, with point\_a being the starting point of the bone tunnel.

**input-point\_b\_index**

Type: Int

Description:

Points are selected using the index, with point\_b being the endpoint of the bone tunnel.

**Output-point\_a**

Type: Point

Description:


point\_a is the starting point of the bone tunnel.

**Output-point\_b**

Type: Point

Description:

point\_b is the endpoint of the bone tunnel.



point_a	tunnel_make	danger_area_of_a
point_b		danger_area_of_b
danger_area		tunnel_line
tunnel_length		cylinder_tunnel
tunnel_width		

## Tunnel make

The input points and set values are used to generate the danger zones and bone tunnels at each point.

### **Input-point\_a**

Type: Point

Description:

point\_a is the starting point of the bone tunnel.

### **Input-point\_b**

Type: Point

Description:

point\_b is the endpoint of the bone tunnel.

### **Input-danger\_area**

Type: Float

Description:

The input value is the radius length of the danger zone.

### **Input-tunnel\_length**

Type: Float

Description:

The input is the length of the bone marrow canal.

### **Input-tunnel\_width**

Type: Float

Description:

The input is the width of the bone marrow canal.

### **Output-danger\_area\_of\_a**

Type: Geometry

Description:

The danger zone generated based on a as the center and the area radius.

### **Output-danger\_area\_of\_b**

Type: Geometry

Description:

The danger zone generated based on b as the center and the area radius.

### Output-tunnel\_line

Type: Line

Description:


The bone tunnel guide line generated based on the point locations inside the bone.

### Output-cylinder\_tunnel



Type: Geometry

Description:

The bone tunnel guide line generated based on the point locations and extension length outside the bone.



**\*All of these outputs should be connected with synchronized geometry in fologram, in an attempt to facilitate AR display and synchronization with HoloLens.**



### distance get

This is used to calculate the distance between the bone drill bit point and the nearest target point among multiple bone tunnel target points.

### Input-spindle\_start\_point

Type: Point

Description:

This point is the drill bit point of the bone drill.

### Input-target\_points

Type: Point

Description:


This point is the target point, which can be a single point or a collection of multiple points.

### Output-distance\_value



Type: Value

Description:

The output value is the distance between the nearest target point and the bone drill bit point.



\*This output should be connected with data visualization.



## Data Visualization

Convert the calculated value into a visualizable text\_entity.

### Input-value

Type: value

Description:

The input value is the calculated distance or angle.



### Input-position

Type: Point

Description:

The location where the text\_entity needs to be generated.

### Input-height

Type: Value

Description:

The size of the text\_entity to be generated.

### Input-angle

Type: Value(0-360)

Description:


Rotate the text\_entity around the generated point location as the center.

### Output-text\_entity



Type: Entity

Description:

Visualize the value as an entity.



\*All of these outputs should be connected with synchronized geometry in fologram, in an attempt to facilitate AR display and synchronization with HoloLens.



degree\_get

By picking the start and end points of the drill bit and the start and end points of the bone tunnel, calculate the angle between the two straight lines of the bone drill and the bone tunnel to facilitate calibration.

#### Input-spindle\_start\_point

Type: Point

Description:

This point is the drill bit point of the bone drill.

#### Input-spindle\_end\_point

Type: Point

Description:

This point is the drill end point of the bone drill.

#### Input-point\_a

Type: Point

Description:

point\_a is the starting point of the bone tunnel.

#### Input-point\_b

Type: Point

Description:


point\_b is the endpoint of the bone tunnel.

#### Output-degree



Type: value

Description:

Calculate the angle between the two straight lines of the bone drill and the bone tunnel.



\*This output should be connected with data visualization.



### surgery line path

This component is used to generate the shortest path for the bone drill to reach the nearest target point.

#### **Input-spindle\_start\_point**

Type: Point

Description:

This point is the drill bit point of the bone drill.

#### **Input-target\_points**

Type: Point

Description:


This point is the target point, which can be a single point or a collection of multiple points.

#### **Output-line**



Type: line

Description:

Generate the shortest path to the nearest target point.



**\*This output should be connected with synchronized geometry**



### surgery curve path

This component is used to generate a curved path for the bone drill to reach the nearest target point, which helps the user calibrate the angle and path.

#### **Input-spindle\_start\_point**

Type: Point

Description:

This point is the drill bit point of the bone drill.

#### **Input-spindle\_end\_point**

Type: Point

Description:

This point is the drill end point of the bone drill.

#### **Input-target\_points**

Type: Point

Description:


This point is the target point, which can be a single point or a collection of multiple points.

#### **Output-curve**


Type: curve

Description:

Generate a Bézier curve to the nearest target point.



**\*This output should be connected with synchronized geometry**



### Extract the bone drill point location

This is a cluster file used to extract the preset start point and end point of the bone drill. We will include it along with the bone drill brep file in the appendix to ensure proper usage.

#### **Input-Brep**

Type: Brep

Type:

Description:

Preset bone drill.

#### **Output-spindle\_start\_point**

Type: Point

Description:

Output the bone drill bit point location.

#### **Output-spindle\_end\_point**

Type: Point


Description:

Output the bone drill end point location.

---

## Usage ExamplesDetailed guide of functionality.

<https://youtu.be/z2X9pDw00ao>



---

## Results From Using Our Tool

We have achieved the precise localization of bone tunnels, which significantly enhances recovery from joint diseases and improves the success rate of surgeries. In this context, MR-based precise bone puncture surgery navigation plays a crucial role in assisting surgical procedures. Unlike traditional methods, we have leveraged the capabilities of MR-based surgical navigation, enabling surgeons, with the aid of MR-SNS, to accurately locate the ideal tunnel positions planned preoperatively and precisely drill the bone tunnels.

Additionally, we have explored the potential of generating an optimal path in static bone puncture surgery, which can notably reduce the cognitive load on surgeons during prolonged operations. This provides a key auxiliary tool within AR's precise navigation system, supporting orthopedic surgeons in performing accurate surgeries.


---

## limitation


Additionally, we have explored the potential of generating an optimal path in static bone puncture surgery, which can notably reduce the cognitive load on surgeons during prolonged operations. This provides a key auxiliary tool within AR's precise navigation system, supporting orthopedic surgeons in performing accurate surgeries.

---

## Bugs



- (1) **bone\_to\_points:** This component runs well in the original Python code, but after compilation, the output point type is incorrect. We have tried multiple versions of the code, but it still does not work.



- (2) **tunnel point select:** This component runs well in the original Python code, but after compilation, it fails to output points. We have tried multiple times, but the code still does not run correctly.